

Product Data Sheet


Product Name *3-D Life* Thioglycerol

Catalog Number T10-3

Description Thioglycerol (CAS-No. 96-27-5) can be covalently coupled with its thiol group to SH-reactive groups of Dextran or PVA polymers to modify *3-D Life* Hydrogels. For more information and instructions, please consult the *3-D Life* Hydrogels User Guide on www.cellendes.com.

Quantity Allows modification of e.g. 9 ml *3-D Life* Hydrogel with 0.4 mmol/L thioglycerol.

Components

Material	Quantity	Concentration of reactive groups	Storage
 Thioglycerol	180 µl	20 mmol/L	-20°C or lower

All materials are filter-sterilized.

Intended for research use only. Not for use in human therapeutic or diagnostic applications.