

Product Data Sheet

Product Name 10x CB (pH 7.2) Phenol Red-free

Catalog Number B21-3

Description *3-D Life* 10x CB-Buffer (pH 7.2) Phenol Red-free is a component of the *3-D Life* Hydrogel system. This buffer does not contain Phenol Red to avoid interferences with light microscopy. It is a 10x concentrated salt and buffer solution which is diluted to physiological strength when used in the *3-D Life* Hydrogel system. It is used for the formation of slow gelling *3-D Life* Hydrogels where it provides a cell compatible environment at pH 7.2.

For instructions, please consult the General Protocol GP-2 "Preparation of *3-D Life* Slow Gelling Hydrogels" and the *3-D Life* Hydrogels User Guide on www.cellendes.com.

Quantity Sufficient for at least 6 ml *3-D Life* Hydrogel

Components

Material	Quantity	Concentration of reactive groups	Storage
 10x CB (pH 7.2) Phenol Red-free	600 µl	n/a	Short term (≤2 months): 4°C Long term: -20°C to -80°C

All materials are filter-sterilized.

Intended for research use only. Not for use in human therapeutic or diagnostic applications.